

Universal Primer List

Primer Name	Sequence(5'-3')	Name	Sequence(5'-3')
1492R	TACGGYTACCTTGTTACGACTT	M13R-pUC	CAGGAAACAGCTATGAC
27F	AGAGTTTGATCMTGGCTCAG	MT_Forward	CATCTCAGTGCAACTAAA
35S-A	AAGGGTCTTGCGAAGGATAG	pBacPAC-RP	GTCTGTAAATCAACAACGC
35S-B	AGTGGAAAAGGAAGGTGGCT	pBAD-F	ATGCCATAGCATTTTTATCCA
518F	CCAGCAGCCGCGTAATACG	pBAD-FP	ATGCCATAGCATTTTTATCC
800R	TACCAGGGTATCTAATCC	pBAD-R	GATTTAATCTGTATCAGG
AD_Reverse	AGATGGTGCACGATGCACAG	pDONOR-FP	TAACGCTAGCATGGATCTC
a-Factor	TACTATTGCCAGCATTGCTGC	pEGFP_N	CCGTCCAGCTCGACCAG
AOX1_Forward	GACTGGTTCCAATTGACAAGC	pEGFP-FP	TTTAGTGAACCGTCAGATC
AOX1_Reverse	GCAAATGGCATTCTGACATCC	pEGFP-RP	AACAGCTCCTCGCCCTTG
BGH-R	TAGAAGGCACAGTCGAGG	pESP-RP	TCCAAAAGAAGTCGAGTGG
Bluescript_KS	TCGAGGTCGACGGTATC	pET-24a	GGGTTATGCTAGTTATTGCTCAG
Bluescript_SK	CGCTCTAGAAGTAGTGGATC	pET-RP	CTAGTTATTGCTCAGCGG
CMV-F	CGCAAATGGGCGGTAGGCGTG	pFastBac_Forward	GGATTATTCATACCGTCCCA
CYC1_Reverse	GCGTGAATGTAAGCGTGAC	pFastBac_Reverse	CAAATGTGGTATGGCTGATT
DsRed1-C	AGCTGGACATCACCTCCCACAACG	pGEX3	GGAGCTGCATGTGTCAGAGG
DsRed1-N	GTAAGGAACTGGGGGGACAG	pGEX5	GGCAAGCCACGTTTGGTG
EBV-RP	GTGGTTTGTCCAAACTCATC	pJET1_2F	CGACTCACTATAGGGAGAGCGGC
EGFP-C	CATGGTCTGCTGGAGTTCGTG	pJET1_2R	AAGAACATCGATTTTCCATGGCAG
EGFP-CF	AGCACCCAGTCCGCCCTGAGC	pMalE	TCAGACTGTGATGAAGC
EGFP-CR	CGTCCATGCCGAGAGTG	pQE-F	CCCGAAAAGTGCCACCTG

EGFP-N	CGTCGCCGTCCAGCTCGACCAG	pQE-R	GTTCTGAGGTCATTACTIONG
EGFP-NR	CGTCGCCGTCCAGCTC	pREP-fwd_primer	GCTCGATACAATAAACGCC
GAL1_Forward	AATATACCTCTACTTTAACGTC	pRH_Forward	CTGTCTCTATACTCCCCTATAG
Gal4AD	TACCACTACAATGGATG	pRH_Reverse	CAAAATCAATAGTTACTATCGC
GLprimer1	TGTATCTTATGGTACTGTAACG	pTrcHis_Forward	GAGGTATATATTAATGTATCG
GLprimer2	CTTTATGTTTTGGCGTCTTCCA	QE_Promoter	CCGAAAAGTGCCACCTG
HCO2198	TAAACTTCAGGGTGACCAAAAAATCA	RVprimer3	CTAGCAAAATAGGCTGTCCC
ITS1	TCCGTAGGTGAACCTGCGG	RVprimer4	GACGATAGTCATGCCCCGCG
ITS2	GCTGCGTTCTTCATCGATGC	SP6	ATTTAGGTGACACTATAG
ITS3	GCATCGATGAAGAACGCAGC	STag_18mer_Primer	GAACGCCAGCACATGGAC
ITS4	TCCTCCGCTTATTGATATGC	SV40-pArev	CCTCTACAAATGTGGTATGG
ITS5	GGAAGTAAAAGTCGTAACAAGG	SV40-Promoter	GCCCCTAACTCCGCCCATCC
KAN2-FP	ACCTACAACAAAGCTCTCATCAACC	T3	ATTAACCCTCACTAAAG
KAN2-RP	GCAATGTAACATCAGAGATTTTGAG	T7	AATACGACTCACTATAG
LCO1490	GGTCAACAAATCATAAAGATATTGG	T7_EEV	ATGTCGTAATAACCCCGCCCCG
M13F	GTA AACGACG GCCAGT	T7promoter	TAATACGACTCACTATAGGG
M13-FP	TGTA AACGACG GCCAGT	T7terminator	GCTAGTTATTGCTCAGCGG
M13F-pUC	GTTTTCCAGTCACGAC	U-19mer_Primer	GTTTTCCAGTCACGACGT
M13R	GCGGATAACAATTTACACAGG		